



Asahipak ODP2 HP

Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel:041 748 50 60; info@infochroma.ch

### Polymer Column Asahipak ODP2-HP

- for the separation of high polar substances
- good selectivity and retention for alkaline and hydrophile substances
- pH 3 - 12

Phase	Plate Number	Particle Size	ID x Length	Column Product Nr.	Price/Pack	Vorsäulen Bestell Nr.	Price/Pack
Asahipak ODP2 HP-4B	3.500	5 $\mu$ m	4.6 x 50 mm	8920-ASpak-ODP2-HP-4B	790.00		
Asahipak ODP2 HP-4D	11.000	5 $\mu$ m	4.6 X 150 mm	8920-ASpak-ODP2-HP-4D	880.00		
Asahipak ODP2 HP-4E	17.000	5 $\mu$ m	4.6 x 250 mm	8920-ASpak-ODP2-HP-4E	910.00		
Asahipak ODP2 HPG-4A		5 $\mu$ m	4.6 x 10 mm			8920-ASpak-ODP2-HPG-4A	550.00
Asahipak ODP2 HP-2B	2.500	5 $\mu$ m	2.0 x 50 mm	8920-ASpak-OPD2-HP-2B	790.00		
Asahipak ODP2 HP-2D	9.000	5 $\mu$ m	2.0 x 150 mm	8920-ASPak-ODP2-HP-2D	880.00		
Asahipak ODP2 HPG-2A		5 $\mu$ m	2.0 x 10 mm			8920-ASpak-ODP2-HPG-2A	550.00

all prices net, excl. VAT, excl. postage and packaging (CHF 12.00); prices subject to change